



	Material	Color(s)	(BSR) Profile	Absorption Rate
VISORB®	Braided	Violet or Undyed	65% @ 14 days	60-90 days








RAZOR POINT™

PRECISION COSMETIC CONVENTIONAL CUTTING NEEDLES 

Needle	Material	Color	Length	Suture Size						
				7-0	6-0	5-0	4-0	3-0	2-0	0
NPC-1 13mm 3/8 circle 	Braided	Undyed	18" (45cm)			834A				
NPC-3 16mm 3/8 circle 	Braided	Undyed	18" (45cm)				845A			

RAZOR POINT™

PRECISION REVERSE CUTTING NEEDLES 

Needle	Material	Color	Length	Suture Size						
				7-0	6-0	5-0	4-0	3-0	2-0	0
NP-1 11mm 3/8 circle 	Braided	Undyed	18" (45cm)		489A*	490A				
NP-3 13mm 3/8 circle 	Braided	Undyed	18" (45cm)		492A	493A	494A			
		Violet	18" (45cm)	463A*		464A				
NPS-3 16mm 3/8 circle 	Braided	Undyed	18" (45cm)			500A				
NPS-2 19mm 3/8 circle 	Braided	Undyed	18" (45cm)			495A	496A	497A		
			30" (75cm)				426A	427A	428A	
		Violet	18" (45cm)				513A*			
NPS-1 24mm 3/8 circle 	Braided	Undyed	18" (45cm)					683A		
			30" (70cm)			935A	936A*			
NP-2 8mm 1/2 circle 	Braided	Undyed	18" (45cm)			503A	504A			
NPS-4 16mm 1/2 circle 	Braided	Undyed	18" (45cm)			506A	507A			

* Call for availability

	Material	Color(s)	(BSR) Profile	Absorption Rate
VISORB®	Braided	Violet or Undyed	65% @ 14 days	60-90 days

REVERSE CUTTING NEEDLES ▾

Needle	Material	Color	Length	Suture Size							
				5-0	4-0	3-0	2-0	0	1	2	
NC-3 13mm 3/8 circle	Braided	Undyed	18" (45cm)	U385A	U386A*						
		Violet	18" (45cm)	385A	386A						
NFS-2 19mm 3/8 circle	Braided	Violet	18" (45cm)	391A	392A	393A					
		Undyed	30" (75cm)	421A	422A	423A					
		Violet	30" (75cm)		397A	398A					
NFS-1 24mm 3/8 circle	Braided	Undyed	30" (75cm)				443A*				
		Violet	30" (75cm)		451A	452A	453A				
NFSL 30mm 3/8 circle	Braided	Undyed	30" (75cm)				589A				
		Violet	30" (75cm)				586A	587A			
NX-1 23mm 1/2 circle	Braided	Violet	30" (75cm)			460A	461A				
NCP-1 37mm 1/2 circle	Braided	Undyed	30" (75cm)				266A				
		Violet	30" (75cm)				466A	467A			



TAPER POINT NEEDLES ○

Needle	Material	Color	Length	Suture Size							
				5-0	4-0	3-0	2-0	0	1	2	
NCT-3 22mm 1/2 circle	Braided	Violet	30" (75cm)			327A	328A				
NSH-1 22mm 1/2 circle	Braided	Violet	30" (75cm)		310A	311A					
NCT-2 26mm 1/2 circle	Braided	Undyed	30" (75cm)				269A				
		Violet	30" (75cm)			332A	333A	334A			
NSH 26mm 1/2 circle	Braided	Undyed	18" (45cm)			416A	417A				
		Violet	30" (75cm)		315A	316A	317A				



* Call for availability

	Material	Color(s)	(BSR) Profile	Absorption Rate
VISORB®	Braided	Violet or Undyed	65% @ 14 days	60-90 days

TAPER POINT NEEDLES

Needle	Material	Color	Length	Suture Size						
				5-0	4-0	3-0	2-0	0	1	2
NCT-1 37mm 1/2 circle 	Braided	Undyed	30" (75cm)				259A*	260A*		
		Violet	30" (75cm)			338A	339A	340A	341A	
		Undyed	36" (90cm)				945A*	946A*		
		Violet	36" (90cm)				345A	346A*		
NCT 40mm 1/2 circle 	Braided	Undyed	30" (75cm)					280A*	281A*	282A
		Violet	30" (75cm)				351A*	352A*	353A*	354A

SOFTCUT NEEDLES

Needle	Material	Color	Length	Suture Size						
				7-0	6-0	5-0	4-0	3-0	2-0	0
NV-7 25mm 1/2 circle 	Braided	Violet	30" (75cm)					180A*	181A*	
NV-34 35mm 1/2 circle 	Braided	Violet	30" (75cm)							364A

MULTI-TY PRE-CUT SUTURES

Needle	Material	Color	Length	Suture Size						
				4-0	3-0	2-0	0	1	2	3
NONE	Braided	Undyed	12-18" (45cm)		110A					
		Violet	12-18" (45cm)		904A*	905A*				

* Call for availability